


 **Water Lens**TM

RAPID COVID-19 DETECTION

FAST, ACCURATE & SIMPLE SALIVA TEST

30 minutes to answer how
you can protect your
community



LOW LOD



ACCURATE
RESULTS



REDUCED FALSE
POSITIVES

4265 San Felipe Suite 1100
Houston, Texas 77027
(844) 987-5367
+1 (713) 960-6620
info@waterlensusa.com
www.waterlensusa.com

For research use only.
Not for for use in diagnostic procedures.



Collect Saliva Sample

Use simple buffered mouth-rinse and have patient spit into collection beaker




Load Sample

Wait 5 minutes and load 7 uL of sample into PCR tube



Get Results

Results in as little as 30 minutes in quantitative or low-cost yes/no response



Point of Care Diagnostics

Collect sample and get results as soon as 30 minutes in the field.



Accurate Results

Sensitive specifically to the SARS-2 virus and not sensitive to the RNA of related viruses (MERS) even at greater than 100,000 copies



Low LOD

Detection limit as low as 10 gene copies and proven effective on intact viruses, capable of identifying less than 300 viral particles.

THE FACTS

- Point of care diagnostic; 3 simple steps, **results in 30 minutes!**
- Utilizes **CDC approved primer set** for optimum specificity and sensitivity.
- Operates with saliva samples. **NO nasal swab required!**
- Choose between **quantitative or low cost qualitative** instruments.
- **Proven reliability.** Of 26 known positive samples, 100% tested positive using this method in lab trials.